


# BioMend® & BioMend® Extend™



30 mm x 40 mm

20 mm x 30 mm

15 mm x 20 mm



## ABSORBABLE COLLAGEN MEMBRANES

### 1 Absorbable<sup>1</sup>

Eliminates second-stage surgery for membrane removal, designed to reduce wound trauma and surgical chair time



BioMend Membrane being prepared for placement in the sinus cavity

### 2 Cell - Occlusive

Serves as barrier to prevent epithelial cell migration and allow passage of essential nutrients<sup>1</sup>

### 3 Space - Maintaining

A rigid enough scaffold for tissue regeneration in GTR and GBR procedures



**ZIMMER BIOMET**  
Your progress. Our promise.®

# BioMend® & BioMend® Extend™

THE ABSORBABLE MEMBRANE RIGID ENOUGH TO CREATE AND MAINTAIN SPACE

## CONVENIENT HANDLING

- Tear resistant, suturable, pliable
- Easy to handle even when hydrated
- Rigid enough to create and maintain space

## VERSATILE BARRIER FUNCTION

- Indicated for use in guided tissue and guided bone regeneration procedures
- BioMend Membrane is designed to be absorbed in approximately 8 weeks<sup>1</sup>
- BioMend Extend Membrane is designed to be absorbed in approximately 18 weeks<sup>1</sup>



Exposed buccal wall defect on #3.



Membrane being prepared for placement.



Membrane cut for placement.



Membrane placed in socket.



Sutures placed to close site.



Final aesthetic result.

Clinical photos ©2013 Gary T. Jones, DMD, Sandhills Oral and Maxillofacial, Fayetteville, NC, USA. All rights reserved. Individual results may vary.

## Ordering Information:

| Catalog Number | Description                  | Catalog Number | Description                         |
|----------------|------------------------------|----------------|-------------------------------------|
| 0103Z          | BioMend Membrane, 15 x 20 mm | 0140Z          | BioMend Extend Membrane, 15 x 20 mm |
| 0105Z          | BioMend Membrane, 20 x 30 mm | 0141Z          | BioMend Extend Membrane, 20 x 30 mm |
| 0107Z          | BioMend Membrane, 30 x 40 mm | 0142Z          | BioMend Extend Membrane, 30 x 40 mm |

<sup>1</sup> Data on file with Collagen Matrix, Inc.

Contact us at 1-800-342-5454 or visit [zimmerbiometdental.com](http://zimmerbiometdental.com)

Zimmer Biomet Dental  
Global Headquarters  
4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
Phone: +1-561-776-6700  
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc., and distributed and marketed by Zimmer Biomet Dental (and, in the case of distribution and marketing, its authorized marketing partners). In certain countries, BioMend and BioMend Extend Membranes are manufactured by Collagen Matrix, Inc. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0056 REV A 08/17 ©2017 Zimmer Biomet. All rights reserved.

